



BIO-MEDICAL EQUIPMENT INSPECTIONS

with Asset Management



Out-Patient/Ambulatory Surgical & Medical Centers



WHAT ARE BIO-MED INSPECTIONS?

1

ELECTRICAL SAFETY INSPECTIONS

All cord-connected medical equipment that come in contact with patients or is used in patient vicinity must receive an annual bio-medical electrical safety inspection. This is an inspection in which the current in the ground wire is measured and must comply with UL specifications. The equipment is then stickered with a “passed” or “failed” sticker specifying the date in which it was inspected. After the safety inspection is completed, a report is generated, listing all inspected equipment with the status and measurement.

2

PREVENTIVE MAINTENANCE

Preventive Maintenance or “PM” service is an equipment inspection and functionality test. Due to the vital nature of some medical equipment for extending un-interrupted patient care it is often the best practice to have the equipment periodically serviced to avoid unexpected equipment failures. Such tests may be as simple as a functionality test, whereas other equipment may require more extensive servicing such as, lubrication of moving parts, wear part replacements (gaskets and filters), and/or extensive cleaning of vital areas of the equipment.

3

EQUIPMENT CALIBRATION

Calibration checks and adjustments may be necessary on various types of medical equipment. A calibration service is an adjustment of certain parameters to meet the actual output or meet the manufacturer’s recommended specification of operation. Depending upon the nature of the calibration, such tests often require test equipment.



OUR GUARANTEE:

We guaranty you an inspection that will pass your state accreditation agency inspection. It also is our commitment that our customers are 100% satisfied.

Medicanix Inc.

Customer Facility Name Here



107001

ASSET TRACKING WITH BAR CODES


Medicanix sets the standards in the industry by giving all the equipment that we inspect its own bar code label. This not only saves time and costs on all future inspections, but it also helps our customers keep track of their equipment and relate it to the written report. All equipment information; model, serial number, equipment type, and test data is recorded and given a control number. A label with the control number and its barcode is put on the equipment. This not only reduces all future inspection time immensely, but also provides our customers with a data record of their own equipment.

KNOWLEDGEABLE PERSONNEL

Medicanix technicians are conversant in knowing which equipment needs inspected. You can have peace of mind that we will perform a thorough job that will enable you to pass whatever accreditation inspection you are trying to accomplish. We will not inspect items you do not wish to inspect without prior authorization but will make recommendations of what should and should not be inspected. You are welcome to accompany our techs during the inspection. However this is not necessary as they generally can work on their own.

INSPECTION LABELS

In order to make it easy to see if an item has been missed or not yet inspected, Medicanix changes the inspection label colors every year. These labels are generally a bright color and easy to see. They are however, put in locations that are less of an eye sore. You are welcome to specify the location of your inspection label as long as it can still be seen by a state or accreditation inspector.

		30 Nurney Street Stamford, CT 06902 (800) 937-3724	
Biomedical Safety Inspection			
PASSED		<DATE>	
NEXT		<NEXT DATE>	
ADDITIONAL TECHNICIAN		<TECH>	





MEDICANIX

Medicanix Inc.

30 Nurney St.

Stamford, CT 06902

Phone: 800-937-3724 or 203-324-3711

Fax: 203-487-7423

Email: sales@medicanix.net

Web: www.medicanix.net



ACCREDITATION ASSOCIATION
for AMBULATORY HEALTH CARE, INC.



DOH
STATE OF NEW YORK
DEPARTMENT OF HEALTH